

Food and Drug Administration Office of Regulatory Affairs**Summary Report****For Sample Number:** 764940

13206132

TD Sample Number:**Import Sample Number**

This is an accurate reproduction of the original electronic record as of 12/10/2012

Sample Class: Normal Everyday Sample**Sample Origin:** Domestic**Sample Basis:** Other-Compl**Sample Flag:****Sample Type:** Official**Collecting District:** NOL-DO**Home District:****Orig C/R and Records To:** NOL-DO**Collection PACs:** 56R814**Product Name:** Methylprednisolone Acetate (Glucocorticoid); Human - Rx/Combination Ingredient; Sterile Liquid**Product Description:** Methylprednisolone Acetate 80mg/ml packaged in 1ml clear glass vials. (Cont in Remarks)**Collection Reason:** Collected during investigation of Aspergillus Meningitis outbreak. Analyze for sterility (presence of Aspergillus) and chemical composition.

Lab: NRL	Split Num: 0	Date Received: 10/05/2012	Date Out of Lab: 12/10/2012
District Conclusion:		District Conclusion Made By:	District
Disposition Reason:		Disposition Authorized By:	Disposition Authorized Date:

Performing Org	PAC	LID	PAF	Compliance No	Lab Class-Description	Laboratory Status
NRL-MSB-K	56R814	M	NAR		3 - Adverse Findings	Completed

Lab Conclusion

Viable microorganisms detected in forty eight (48) of seventy eight (78) units examined.

Yeast-like organisms detected in units # 3, 4, 6, 7, 8, 10, 13, 28E, 29A, 29D, 30A, 30D, 30E, 31B, 31D, 32B, 33A, 33B, 34C, 35A, 35E, 36D, 37B, 38B, 38C, 38D, 38E and 40E.

Yeast-like and mold-like organisms detected in units # 12, 28A, 28B, 28D, 29C, 30C, 31C, 32A, 32C, 32D, 33C, 35C, 35D, 37C, 37D, 37E, 39E, 40B and 40C.

Yeast-like and fine cotton candy like, mold-like, non-septated hyphae organisms detected in unit # 36B.

No viable microorganisms detected in thirty (30) of seventy eight (78) units tested (units # 1, 2, 5, 9, 11, 28C, 29B, 29E, 30B, 31A, 31E, 32E, 33D, 33E, 34A, 34B, 34D, 34E, 35B, 36A, 36C, 36E, 37A, 38A, 39A, 39B, 39C, 39D, 40A, 40D).

No viable microorganisms detected from swabs taken from the exterior portion of the rubber stopper of each of seventy eight (78) units cultured in Sab-Dextrose (Sab-Dex) medium.


No obvious foreign material observed under direct visual examination of the seventy eight (78) of seventy eight (78) units individually examined.

All media, equipment, environmental and system controls were negative for growth.

Lab Conclusion Date**Lab Conclusion Made By**

12/10/2012

Braut Taormina, Jasna



Date: 12/10/2012

Page: 1 of 1

FDA_E00277930

Flag:

ANALYST WORKSHEET		1. PRODUCT Methylpred. AC(PF) 80MG/ML Injectable		2. SAMPLE NO. 764940	
3. SEALS ☒ INTACT ☐ NONE ☐ BROKEN		4. DATE RECEIVED 10/5/2012		5. RECEIVED FROM Elvis Leslie	
6. DISTRICT OR LABORATORY NRL					
7. DESCRIPTION OF SAMPLE One brown paper bag officially sealed as "764940, 10/2/2012, Marshall H. Kinsey, Investigator" containing a total of twenty seven (27) intact individuals units (vials), and forty (42) intact foil packs containing 5 vials each of product; for a total of 69 subs which is a total of 237 units. Each vial and foil packs of product were placed in a whirl-pak bag individually. Whirl-pak bags were identified as "764940, 10/02/12 MHR Subs" 1 to 69. Form FDA 525 attached. Sample received at ambient temperature.					
8. NET CONTENTS		NOT APPLICABLE DECLARE/UNITS <u>1 ml</u> ☒ NOT DETERMINED AMOUNT FOUND _____ UNITS EXAMINED % OF DECLARED _____		9. LABELING	
				1 set & 2 sets copies of ORIGINAL SUBMITTED PHOTOCOPIES SUBMITTED NONE	
10. SUMMARY OF ANALYSIS					
A. CONTAINER: Primary container: Clear glass vial with a gray rubber stopper under a crimp metal closure and a white pop off plastic dust cover measuring approximately 3.5 cm in height and 1.0 cm in base diameter. Secondary container: Foil bag with ziplock type closure measuring approximately 12cm H x 11.5cm W. B. LABELING: Commercially manufactured pressure sensitive paper label affixed around glass vial. C. CODE: "Lot#08102012@51 Discard after 2/6/2013" printed on paper label affixed around glass vial. D. PRODUCT: White opaque free flowing suspension. E. ANALYSIS: Examination for viable microorganism. Sterility analysis. Direct visual examination. Exterior rubber stopper for objectionable microorganism. F. METHOD: USP 35-NF 30, 2012 Section <71> Sterility Test. Vol. 1 pp 69-74 Sterility Analytical Manual, 3rd Edition, August 1997, pp. 5-17 See p. 4 for methodology. G. UNITS TESTED: Seventy eight (78) units tested. H. RESULTS: Viable microorganisms detected in forty eight (48) of seventy eight (78) units tested. Yeast-like organisms detected in units # 3, 4, 6, 7, 8, 10, 13, 28E, 29A, 29D, 30A, 30D, 30E, 31B, 31D, 32B, 33A, 33B, 34C, 35A, 35E, 36D, 37B, 38B, 38C, 38D, 38E, 40E. Yeast-like and mold-like organisms detected in units #12, 28A, 28B, 28D, 29C, 30C, 31C, 32A, 32C, 32D, 33C, 35C, 35D, 37C, 37D, 37E, 39E, 40B, 40C. Yeast-like and fine cotton candy like, mold-like, non-septated hyphae organisms detected in unit 36B. No viable microorganisms detected in thirty (30) of seventy eight (78) units tested (Units # 1, 2, 5, 9, 11, 28C, 29B, 29E, 30B, 31A, 31E, 32E, 33D, 33E, 34A, 34B, 34D, 34E, 35B, 36A, 36C, 36E, 37A, 38A, 39A, 39B, 39C, 39D, 40A, 40D). No viable microorganisms detected from swabs taken from the exterior portion of the rubber stopper of each of seventy eight (78) units cultured in Sab-Dextrose (Sab-Dex) medium. No obvious foreign material observed under direct visual examination of the seventyeight (78) of seventyeight (78) units individually examined. All media, equipment, environmental and system controls were negative for growth.					
11. RESERVE SAMPLE: See page # 2					
12. ANALYST SIGNATURE Haydee B. Romero (x) Broke Seal 10/5/12 <i>Haydee B. Romero</i> <i>Jessica Brand - Thormina</i> <i>Daniel J. O'Reilly</i>			13. WORK SHEET CHECK a. BY <i>[Signature]</i> Ph.D. b. DATE 12/10/12		
			14. DATE REPORTED 12/10/12		

GENERAL CONTINUATION SHEET	PRODUCT: Methylpred. AC(PF) 80MG/ML Injectable	SAMPLE NO. 764940
11. RESERVE SAMPLE: Brown paper bag officially sealed "764940, 10/5/2012, Haydee B. Romero" containing a total of fourteen (14) intact (individuals vials) subs, twenty nine (29) intact (foil packs) subs for a total of 145 units, which represent a total of 43 intact subs which is a total of 159 intact units. There is also a total of thirteen (13) empty (individuals vials) subs, thirteen (13) open (foil packs) subs for a total of 65 empty units; which represent a total of 26 open empty subs which is a total of 78 units of which the entire content has been used in analysis. Each unit is identified as "764940 10/5/12 HR". Two (2) whirl pak bags officially sealed "764940, 10/10/12, Haydee B. Romero, Analyst". Each bag contains eight isolates on Sab-Dextrose agar slant; slants identified as "Sample # 764940 Sub#12-1 (pink colony) 10/10/12 HR", "Sample # 764940 Sub #12-2 (gray colony) 10/10/12 HR", "Sample # 764940 Sub # 29A (pink colony) 10/10/12 HR", "Sample # 764940 Sub # 36B 10/10/12 HR", "Sample # 764940 Sub#37D-1 (pink colony) 10/10/12 HR", Sub #37D-2 (gray colony) 10/10/12 HR", "Sample # 764940 Sub#39E-1 (pink colony) 10/10/12 HR", "Sub #39E-2 (gray colony) 10/10/12 HR". One bag was shipped to CDC on 10/11/12, for further identification and the other bag held within a secured refrigerator at NRL. One label and one empty pack were sacrificed to include with analytical worksheet. Sample kept under lock and key until returned to sample custodian.		
ANALYST(S) <i>Haydee B Romero</i>		PAGE 2 OF 3 PAGES

GENERAL CONTINUATION SHEET	PRODUCT <u>Methylpred. AC (PF) 80 MG / ML</u> <u>Injectable</u>	SAMPLE NO. <u>764940</u>			
STERILITY MEDIA & EQUIPMENT CONTROL SHEET					
Incubation Started: <u>10/5/12</u> By: <u>HZ</u>					
CONTROLS: Environmental Monitoring: <u>Two Blood agar</u> plates exposed for the duration of the analysis. (Maximum time of exposure is 1 hour) Agar Plate Growth: <u>0</u> bacterial colonies <u>0</u> mold colonies <u>0</u> yeast colonies					
Media and Equipment Quality Assurance Numbers and Control Results:					
ITEM	QA #	Results	ITEM	QA #	Results
Thioglycollate Broth	276-12	<u>Neg</u>	Millipore Forceps	N/A	N/A
Soybean Casein Digest Medium	177-12	<u>Neg</u>	Manifolds	N/A	N/A
Rinse Fluid A	N/A	N/A	Syringes R'd 8/25/98 <u>B.D. ICC 256 5/8</u>	Lot 0K345	<u>Neg</u>
Rinse Fluid D	N/A	N/A	Beakers	131-12	N/A
Malt Extract Agar (MEA)	276-12	<u>Neg</u>	Chlortetracycline	242-12	N/A
3% Hydrogen Peroxide R'd 3/27/12	Lot 11524	N/A	Swabs R'd 9/14/11	Lot 2828	<u>Neg</u>
70% Filter Sterilized alcohol	264-11	N/A	Membranes .22microns	N/A	N/A
			Membranes .45microns	N/A	N/A
Wipes	176-12	N/A	10" Forceps	047-12	<u>Neg</u>
10 ml Sabdex Broth	276-12	<u>Neg</u>	Phosphate Buffer	276-12	<u>Neg</u>
System Control	011-11	<u>Neg</u>	Other/Blood Agar BBL	Lot 2220296	<u>Neg</u>
22.5+/-2.5C Incubator <u>1776</u> 32.5+/-2.5 Incubator <u>1958</u> Neg = Negative					
N/A=Not Applicable					
ANALYST(S): <u>Haydee B Romero</u>				Page <u>3</u> of <u>31</u> Pages	

General continuation sheet	Product <u>Methylpred. AC (PF)</u> <u>80N6/ML Insectable</u>	Sample No. <u>764940</u>
Methodology-Direct inoculation Sample preparation , swabbing under dust cap and sterility testing (<u>Vials</u>)		
1. Sample Preparation: The number of unit (s) <u>(Seventy eight (78))</u> to be examined was/were removed from the sample container under the laminar flow hood and all units were placed in sterile bag. The bag was disinfected using 3% H2O2 disinfectant solution. The disinfected bag then placed inside the disinfected stainless steel pass thru of the clean room and left for an hour. Analyst wore sterile latex rubber gloves, sterile sleeve, hair net and a mask during this procedure (M).		
2. Preparation of the analytical area and personnel: Within the gowning area the analyst (s) dressed in a garment (s) consisting of hair net, mask, shoe cover, sterile hood, sterile jumpsuit, sterile booties and sterile latex rubber gloves. All transfers of product to growth medium were aseptically performed within a horizontal laminar flow hood (LFH) and the Bio safety cabinet (BSC) located within the bio-clean room. All surfaces within the LFH and BSC were thoroughly disinfected with <i>Filter sterilized 70% alcohol</i> prior to analysis of the sample (M).		
3. Analytical procedures: All units were placed within the LFH. For each unit the flip off dust cap was aseptically removed and by using sterile cotton swap moisten in sterile <i>PB Buffer</i>, the exposed rubber stopper area was wiped and the swab was aseptically transferred to a tube containing 10 ml of sterile <i>Sab-Dextrose</i> (Sab-Dex) medium (M). All tested unit (s) was/were then transferred to another LFH hood within the clean room for disinfection (M). Each unit was wiped with 3% H2O2 disinfectant solution and left to dry for up to one (1) hour (M). All disinfected unit(s) was/were transferred to the BSC within the clean room (M). Within the BSC within the Bio clean room. These units were again disinfected with <i>Filter sterilized 70% alcohol</i> and left to air dry under the BSC (M). For each unit using sterile <u>1</u> CC syringe <u>0.5</u> ml of the contents were aseptically delivered into SCD medium and <u>0.5</u> ml delivered into Fl. Thio. medium (M). The remaining droplets in the syringe were then aseptically transferred unto a MEA agar plate (M).		
4. Controls: <u>Media:</u> All tubes of Fl. Thio medium, SCD medium , 10 ml <u>Sab-dex</u> medium , blood agar plates, and <u>MEA</u> agar plates were pre-incubated prior to use. One un-inoculated tube of each medium was incubated with the test tubes as closed controls (M). <u>Swab:</u> One sterile swab was aseptically transferred to a tube containing 10 ml of sterile <u>Sab-dex</u> medium. <u>Equipment:</u> Two (2) <u>10" forceps</u> were each dipped into one tube each of SCD and Fl. Thio. media (M). <u>Syringe:</u> Two (2) syringes were used to withdraw media from each of SCD and Fl. Thio. Media and expelled it back to their respective tubes. <u>System control:</u> To monitor the sterility of the analytical technique, a system control was used. <i>One sterile vial of 0.1% peptone water was aseptically opened to expose the rubber stopper by using a sterile 10" forceps, this unit was tested using the same analytical procedure employed in setting up the sample</i> (M).		
5. Environmental: Two <i>blood</i> agar plates were exposed for the duration of the analysis on the right and left sides of the hood working area for each step of the analysis (M). One tube each of Fl. Thio. medium, SCD medium and SabDex medium were opened and exposed for the duration of the analysis as opened controls (M).		
6. Incubation: Tubes of SCD medium were incubated <u>14</u> days at 22.5 +/- 2.5 °C. Tubes of Fl. Thio. medium were incubated <u>14</u> days at 32.5 +/- 2.5 °C. Tubes of Sab-Dextrose medium were incubated 14 days at 30 °C. All tubes were observed periodically for indication of possible growth (i.e. turbidity) (M). The MEA agar plates were incubated at 22.5 +/- 2.5 °C and were observed for growth for up to 14 days (M).		
7. Subculture of presumptive positive tubes: <u>See p. 22 - 26</u>		
() = Analyst performing procedure		
ANALYST(S) <u>Huydee B. Romero</u>		PAGE <u>4</u> OF <u>31</u> PAGES

DAILY OBSERVATIONS

SAMPLE TEST

Date Incubation Started 10/5/12 HCL

SUB Composite	Thio/SCD DAYS INCUBATED														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
1				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page #1 and pages 22-31
2				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
3				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
4				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
5				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
6				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
7				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
8				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
9				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
10				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
11				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
12				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
13				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
N/A															V
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 H	10/10 H	10/11 H	10/12 H	10/13 W	10/14 W	10/15 H	10/16 H	10/17 H	10/18 H	10/19 H	

Upon receiving the sample a direct visual examination was done. No obvious foreign material observed under direct visual examination of 78 units individually examined.

+ = obvious growth observed

(+) = presumptive positive, unable to identify growth due to the turbidity nature of the product, all Thio/SCD tubes were streaked to MacConkey agar, Soybean Casein Digest Agar and H₂O extract agar after 14 days.

W = weekend H = Holidays

ANALYST(S) Haydee B. Romero

Page 5 of 31 Pages

DAILY OBSERVATIONS

SAMPLE TEST

Date Incubation Started 10/5/12

SUB / Composite	DAYS INCUBATED														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
28 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page # 1 and pages 22-31
28 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
28 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
28 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
28 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
29 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
29 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
29 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
29 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
29 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
30 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
30 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
30 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
30 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
30 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 M	10/10 M	10/11 M	10/12 M	10/13 W	10/14 W	10/15 M	10/16 M	10/17 M	10/18 M	10/19 M	

See page 5 for Direct Visual Examination

+ = obvious growth observed

(+) = presumptive positive, unable to identify growth due to turbidity
 nature of the product, all Thio/SCD Tubes were streaked to MacConkey
 agar, Soybean Casein Digest Agar and Malt extract agar after 14 days.

W = weekend

H = Holidays

ANALYST(S)

Haydee B. Roman

Page 6 of 31 Pages

DAILY OBSERVATIONS

SAMPLE TEST

Date Incubation Started 10/5/12 dbr

SUB Composite	Thio/SCD DAYS INCUBATED														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
31 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page 41 and page 22-31
31 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
31 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
31 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
31 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
32 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
32 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
32 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
32 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
32 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
33 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
33 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
33 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
33 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
33 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 M	10/10 M	10/11 M	10/12 M	10/13 W	10/14 W	10/15 M	10/16 M	10/17 M	10/18 M	10/19 M	

See page 5 for Direct Visual Examination

+ = obvious growth observed
 (+) = presumptive positive, unable to identify growth due to the turbidity nature of the product, all Thio/SCD tubes were streaked to MacConkey agar, Sybian Casein Digest Agar and Malt Extract agar after 14 days.
 W = Weekend H = Holidays

ANALYST(S)

Haylee B. Amew

Page 7 of 31 Pages

DAILY OBSERVATIONS

SAMPLE TEST

Date Incubation Started 10/5/12 *br*

SUB / Composite	Thio/SCD DAYS INCUBATED														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
34 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page 11 and pages 22-31
34 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
34 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
34 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
34 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
35 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page 11 and pages 22-31
35 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
35 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
35 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
35 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
36 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
36 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
36 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
36 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
36 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 W	10/10 W	10/11 W	10/12 W	10/13 W	10/14 W	10/15 W	10/16 W	10/17 W	10/18 W	10/19 W	

See pages 5 for Direct Visual Examination

+ = obvious growth observed

(+) = presumptive positive, unable to identify growth due to the turbid nature of the product, all Thio/SCD tubes were streaked to MacConkey agar, Soybean Casein Digest agar and H₂O₂ T₂ agar after 14 days

W = Weekend

H = Holidays

ANALYST(S)

Haylee B. Rouse

Page 8 of 31 Pages

DAILY OBSERVATIONS

SAMPLE TEST

Date Incubation Started 10/5/12/12

SUB / Composite	Thio/SCD DAYS INCUBATED														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
37 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page #1 and pages 22-31
37 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
37 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
37 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
37 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
38 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
38 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
38 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
38 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
38 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
39 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
39 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
39 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
39 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
39 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
DATE INITIAL	10/6/12 W	10/7/12 W	10/8/12 H	10/9/12 W	10/10/12 W	10/11/12 W	10/12/12 W	10/13/12 W	10/14/12 W	10/15/12 W	10/16/12 W	10/17/12 W	10/18/12 W	10/19/12 W	

See page 5 for Direct Visual Examination

+ = obvious growth observed.

(+) = presumptive positive, unable to identify growth due to the turbidity nature of the inoculum, all Thio/SCD tubes were streaked to MacConkey agar, SBA, and casein digest agar and H₂O₂ test after 14 days

W = Weekend H = Holidays

ANALYST(S) Haydee B Romero

Page 9 of 31 Pages

DAILY OBSERVATIONS

SAMPLE TEST

Date Incubation Started *10/5/12*

SUB Composite	Thio/SCD DAYS INCUBATED														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
40 A				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	See page #1 and pages 22-31 ↓
40 B				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
40 C				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
40 D				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
40 E				(+)	(+)	(+)	(+)			(+)	(+)	(+)	(+)	(+)	
System Control				-	-	-	-			-	-	-	-	-	No Growth ↓
Open Media				-	-	-	-			-	-	-	-	-	
Closed Media				-	-	-	-			-	-	-	-	-	
Syringes				-	-	-	-			-	-	-	-	-	
10" forceps				-	-	-	-			-	-	-	-	-	
1/A															↓
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 W	10/10 W	10/11 H	10/12 H	10/13 W	10/14 W	10/15 W	10/16 W	10/17 W	10/18 W	10/19 W	

See page 5 for Direct Visual Examination

+ = obvious growth observed

(+) = presumptive positive, unable to identify growth due to the turbid nature of the product, all Thio/SCD tubes were streaked to MacConkey agar, trypticase soy agar and H₂O extract

W = Weekend

H = Holiday

agar after 14 days

ANALYST(S)

*Haydee B. Romero*Page *10* of *31* Pages

DAILY OBSERVATIONS

From Swabs Taken from the exterior portion of the rubber stoppers of each sub.

SAMPLE TEST

Date Incubation Started *10/5/12 BL*

SUB / Composite	swab in Sab Dex broth														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
1				-	-	-	-			-	-	-	-	-	No Growth ↓
2				-	-	-	-			-	-	-	-	-	
3				-	-	-	-			-	-	-	-	-	
4				-	-	-	-			-	-	-	-	-	
5				-	-	-	-			-	-	-	-	-	
6				-	-	-	-			-	-	-	-	-	
7				-	-	-	-			-	-	-	-	-	
8				-	-	-	-			-	-	-	-	-	
9				-	-	-	-			-	-	-	-	-	
10				-	-	-	-			-	-	-	-	-	
11				-	-	-	-			-	-	-	-	-	
12				-	-	-	-			-	-	-	-	-	
13				-	-	-	-			-	-	-	-	-	
N/A															
↓															
DATE INITIAL	<i>10/6 W</i>	<i>10/7 BL</i>	<i>10/8 H</i>	<i>10/9 BL</i>	<i>10/10 BL</i>	<i>10/11 BL</i>	<i>10/12 BL</i>	<i>10/13 W</i>	<i>10/14 W</i>	<i>10/15 BL</i>	<i>10/16 BL</i>	<i>10/17 BL</i>	<i>10/18 BL</i>	<i>10/19 BL</i>	

W= weekend L= Leave H= Holidays O= other (Explain)

- = No Growth

ANALYST(S)

*Haydee & Romero*Page *11* of *31* Pages

DAILY OBSERVATIONS

From swabs Taken from the exterior portion of the rubber stopper of each sub.

SAMPLE TEST

Date Incubation Started 10/5/12 OR

(SUB) Composite	swab in Sab Dex broth														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
28 A				-	-	-	-			-	-	-	-	-	No Growth
28 B				-	-	-	-			-	-	-	-	-	
28 C				-	-	-	-			-	-	-	-	-	
28 D				-	-	-	-			-	-	-	-	-	
28 E				-	-	-	-			-	-	-	-	-	
29 A				-	-	-	-			-	-	-	-	-	
29 B				-	-	-	-			-	-	-	-	-	
29 C				-	-	-	-			-	-	-	-	-	
29 D				-	-	-	-			-	-	-	-	-	
29 E				-	-	-	-			-	-	-	-	-	
30 A				-	-	-	-			-	-	-	-	-	
30 B				-	-	-	-			-	-	-	-	-	
30 C				-	-	-	-			-	-	-	-	-	
30 D				-	-	-	-			-	-	-	-	-	
30 E				-	-	-	-			-	-	-	-	-	✓
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 M	10/10 M	10/11 M	10/12 M	10/13 W	10/14 W	10/15 M	10/16 M	10/17 M	10/18 M	10/19 M	

W= weekend

L= Leave

H= Holidays

O= other (Explain)

- = NO Growth

ANALYST(S)

Haydee BrumerPage 12 of 31 Pages

DAILY OBSERVATIONS

From swabs taken from the exterior portion of the rubber stopper of each sub.

SAMPLE TEST

Date Incubation Started 10/5/12 DN

(SUB) / Composite	swab in Sab Dex broth														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
31 A				-	-	-	-			-	-	-	-	-	No Growth
31 B				-	-	-	-			-	-	-	-	-	
31 C				-	-	-	-			-	-	-	-	-	
31 D				-	-	-	-			-	-	-	-	-	
31 E				-	-	-	-			-	-	-	-	-	
32 A				-	-	-	-			-	-	-	-	-	
32 B				-	-	-	-			-	-	-	-	-	
32 C				-	-	-	-			-	-	-	-	-	
32 D				-	-	-	-			-	-	-	-	-	
32 E				-	-	-	-			-	-	-	-	-	
33 A				-	-	-	-			-	-	-	-	-	
33 B				-	-	-	-			-	-	-	-	-	
33 C				-	-	-	-			-	-	-	-	-	
33 D				-	-	-	-			-	-	-	-	-	
33 E				-	-	-	-			-	-	-	-	-	
DATE INITIAL	10/6 W	10/7 W	10/8 H	10/9 DN	10/10 DN	10/11 DN	10/12 DN	10/13 W	10/14 W	10/15 O	10/16 DN	10/17 DN	10/18 DN	10/19 DN	

W= weekend L= Leave H= Holidays O= other (Explain) - = No Growth

ANALYST(S)

Haydee B. RomeroPage 13 of 31 Pages

DAILY OBSERVATIONS

From swabs taken from the *external* portion of the rubber stopper of each sub

SAMPLE/TEST

Date Incubation Started *10/5/12 BL*

SUB Composite	swab in Sab Dex broth														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
34 A				-	-	-	-			-	-	-	-	-	No Growth
34 B				-	-	-	-			-	-	-	-	-	
34 C				-	-	-	-			-	-	-	-	-	
34 D				-	-	-	-			-	-	-	-	-	
34 E				-	-	-	-			-	-	-	-	-	
35 A				-	-	-	-			-	-	-	-	-	
35 B				-	-	-	-			-	-	-	-	-	
35 C				-	-	-	-			-	-	-	-	-	
35 D				-	-	-	-			-	-	-	-	-	
35 E				-	-	-	-			-	-	-	-	-	
36 A				-	-	-	-			-	-	-	-	-	
36 B				-	-	-	-			-	-	-	-	-	
36 C				-	-	-	-			-	-	-	-	-	
36 D				-	-	-	-			-	-	-	-	-	
36 E				-	-	-	-			-	-	-	-	-	
DATE INITIAL	<i>10/6 W</i>	<i>10/7 W</i>	<i>10/8 H</i>	<i>10/9 BL</i>	<i>10/10 BL</i>	<i>10/11 BL</i>	<i>10/12 BL</i>	<i>10/13 W</i>	<i>10/14 W</i>	<i>10/15 BL</i>	<i>10/16 BL</i>	<i>10/17 BL</i>	<i>10/18 BL</i>	<i>10/19 BL</i>	

W= weekend

L= Leave

H= Holidays

O= other (Explain)

-- No Growth

ANALYST(S)

*Haydee B Romero*Page *14* of *31* Pages

DAILY OBSERVATIONS

From swabs Taken from the exterior portion of the rubber stopper of each sub.

SAMPLE TEST

Date Incubation Started *10/5/12*

(SUB) Composite	swab in Sab Dex broth														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
37 A				-	-	-	-			-	-	-	-	-	No Growth
37 B				-	-	-	-			-	-	-	-	-	
37 C				-	-	-	-			-	-	-	-	-	
37 D				-	-	-	-			-	-	-	-	-	
37 E				-	-	-	-			-	-	-	-	-	
38 A				-	-	-	-			-	-	-	-	-	
38 B				-	-	-	-			-	-	-	-	-	
38 C				-	-	-	-			-	-	-	-	-	
38 D				-	-	-	-			-	-	-	-	-	
38 E				-	-	-	-			-	-	-	-	-	
39 A				-	-	-	-			-	-	-	-	-	
39 B				-	-	-	-			-	-	-	-	-	
39 C				-	-	-	-			-	-	-	-	-	
39 D				-	-	-	-			-	-	-	-	-	
39 E				-	-	-	-			-	-	-	-	-	
DATE INITIAL	<i>10/6 W</i>	<i>10/7 W</i>	<i>10/8 H</i>	<i>10/9 M</i>	<i>10/10 M</i>	<i>10/11 M</i>	<i>10/12 M</i>	<i>10/13 W</i>	<i>10/14 W</i>	<i>10/15 M</i>	<i>10/16 M</i>	<i>10/17 M</i>	<i>10/18 M</i>	<i>10/19 M</i>	

W= weekend

L= Leave

H= Holidays

O= other (Explain)

- = No Growth

ANALYST(S)

*Huydee B. Romero*Page *15* of *31* Pages

DAILY OBSERVATIONS

From Swabs taken from the exterior portion of the rubber stoppers of each sub.

SAMPLE TEST

Date Incubation Started *10/5/12*

SUB Composite	swab in Sab Dex broth														Final Results
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
40 A				-	-	-	-			-	-	-	-	-	No Growth
40 B				-	-	-	-			-	-	-	-	-	
40 C				-	-	-	-			-	-	-	-	-	
40 D				-	-	-	-			-	-	-	-	-	
40 E				-	-	-	-			-	-	-	-	-	
Swab control				-	-	-	-			-	-	-	-	-	
Open control				-	-	-	-			-	-	-	-	-	
Closed control				-	-	-	-			-	-	-	-	-	
System Control				-	-	-	-			-	-	-	-	-	
<i>1/A</i>															
DATE INITIAL	<i>10/6 W</i>	<i>10/7 W</i>	<i>10/8 H</i>	<i>10/9 BL</i>	<i>10/10 BL</i>	<i>10/11 BL</i>	<i>10/12 BL</i>	<i>10/13 W</i>	<i>10/14 W</i>	<i>10/15 BL</i>	<i>10/16 BL</i>	<i>10/17 BL</i>	<i>10/18 BL</i>	<i>10/19 BL</i>	

W= weekend

L= Leave

H= Holidays

O= other (Explain)

- = No Growth

ANALYST(S)

*Haydee B. Rouen*Page *16* of *31* Pages

GENERAL CONTINUATION SHEET		PRODUCT <i>Methy/pred. AC (PF) 80 H6/HL</i> <i>Injectable</i>		SAMPLE NO. <i>764940</i>	
Direct inoculation of product on agar media				N/A=Not Applicable	
SUB #	Agar type used	Plates read <i>10/10/12/HL</i> <i>10/9/12/HL</i>	Gram stain <i>10/9/12/HL</i>	Subculture on <i>10/10/12/HL</i>	Biochemical identification <i>N/A</i>
<i>10/5-12/HL</i>					
1	MEA	—	N/A		
2		—	N/A		
3		+	Y		<i>Yeast like organisms</i>
4		+	* *		
5		—	N/A		
6		+	* *		
7		+	↓		
8		+	↓		
9		—	N/A		
10		+	* *		
11		—	N/A		
12		+	Y & M	* <i>Sab-dex Skirt</i>	<i>Yeast like and Mold like organisms</i>
13		+	* *		
28 A		+	↓		
28 B		+	↓		
28 C		—	N/A		
28 D		+	* *		
28 E	↓	+	↓		
N/A ↓					

ANALYST(S) *Hauydee B. Romero*

PAGE *170B3* / PAGES

Key: Y = Yeast-like cells

M = Septated Hyphae

GENERAL CONTINUATION SHEET		PRODUCT <i>Me Thyl pred. AC (PF) 80MG/HL</i>		SAMPLE NO. <i>764940</i>		
Direct inoculation of product on agar media						
SUB #	Agar type used	Plates read <i>10/9/12 HZ</i> <i>10/9/12 HZ</i>	Gram stain <i>10/9/12 HZ</i>	Subculture on <i>10/10/12 HZ</i>	Biochemical identification <i>N/A</i>	Results
<i>10/5/12 HZ</i> 29 A	MEA	+	Y	<i>Sab-dex Skint</i>		<i>yeast-like organisms</i>
29 B		-	N/A			
29 C		+	* ↓			
29 D		+	↓			
29 E		-	N/A			
30 A		+	* *			
30 B		-	N/A			
30 C		+	* *			
30 D		+	↓			
30 E		+	↓			
31 A		-	N/A			
31 B		+	* *			
31 C		+	↓			
31 D		+	↓			
31 E		-	N/A			
32 A		+	* *			
32 B		+	↓			
32 C		+	↓			
32 D		+	↓			
32 E	↓	-	N/A			

ANALYST(S) *Huydee B Romero*

PAGE *18* OF *1* PAGES

Key: Y = Yeast-like cells

M = Septated Hyphae

GENERAL CONTINUATION SHEET		PRODUCT <i>Methylpred. AC (PF) 804G/HL</i>		SAMPLE NO. 764940		
Direct inoculation of product on agar media					N/A=Not Applicable	
SUB #	Agar type used	Plates read <i>10/9/12 GR</i> <i>10/9/12 GR</i>	Gram stain <i>10/9/12 GR</i>	Subculture on <i>10/10/12 GR</i>	Biochemical identification N/A	Results
33 A	MEA	+	* *			
33 B		+	↓			
33 C		+	↓			
33 D		-	N/A			
33 E		-	N/A			
34 A		-	N/A			
34 B		-	N/A			
34 C		+	* *			
34 D		-	N/A			
34 E		-	N/A			
35 A		+	* *			
35 B		-	N/A			
35 C		+	* *			
35 D		+	↓			
35 E		+	↓			
36 A		-	N/A			
36 B		+	MC	Sub-dex Slant		Mold-like organisms
36 C		-	N/A			
36 D		+	* *			
36 E	✓	-	N/A			

ANALYST(S) <i>Haydee B Rouse</i>	PAGE/9 OF 3 / PAGES
----------------------------------	---------------------

Key: Y = Yeast-like Cells

M = Septated Hyphae
MC = non-septated hyphae

GENERAL CONTINUATION SHEET		PRODUCT <i>Methy</i> pred. AC (PF) 80 MG/ML Injectable		SAMPLE NO. 764940		
Direct inoculation of product on agar media					N/A=Not Applicable	
SUB #	Agar type used	Plates read	Gram stain	Subculture on	Biochemical identification	Results
10/5/12 <i>DL</i>		10/9/12 <i>DL</i> 10/12/12 <i>DL</i>	10/9/12 <i>DL</i>	10/10/12 <i>DL</i>	N/A	
37 A	MEA	-	N/A			
37 B		+	* *			
37 C		+	↓			
37 D		+	Y & M	Sab-dex Slant		yeast-like and mold-like organisms.
37 E		+	* *			
38 A		-	N/A			
38 B		+	* *			
38 C		+	↓			
38 D		+	↓			
38 E		+	↓			
39 A		-	N/A			
39 B		-	N/A			
39 C		-	N/A			
39 D		-	N/A			
39 E		+	Y & M	Sab-dex Slant		yeast-like and mold-like organisms.
40 A		-	N/A			
40 B		+	* *			
40 C		+	↓			
40 D		-	N/A			
40 E	↓	+	* *			

ANALYST(S) *Haylee B Rowen*

PAGE 20 OF 31 PAGES

Key: Y = Yeast-like cells M = Septated Hyphae

[illegible]

PRODUCT Methylpred. AC (PF) 80 MG / ML
Injectable

SAMPLE NO.
764940

Direct inoculation of product on agar media

N/A=Not Applicable

SUB #	Agar type used	Plates read 10/19/12/16/2	Gram stain N/A	Subculture on N/A	Biochemical identification N/A	Results
Blank/Control	MEA	—				No Growth
N/A						
*	A representative portion of positive subs were subcultured to Sab-dex slants, in duplicate sets (Sub 12 = yeast like colony & mold like colony, sub 29A = yeast like colony, sub 36 B = mold-like colony, Sub 37 D = yeast like colony & mold-like colony, Sub 39 E = yeast like colony & mold-like colony) One set for shipment to CDC for further identification and the other set for storage at NRL					
**	Colony morphology same as sub# 3 and sub 29A					
Sab Dex agar slants QA # 279-12						
Microscope 282						

ANALYST(S) Haydee B Romero

PAGE 21 OF 31 PAGES

GENERAL CONTINUATION SHEET				PRODUCT <i>Methylpred. AC (PF) 8046/ML</i>			SAMPLE NO. <i>769940</i>		
Results from Subculture of <u><i>THIO</i></u> Medium on agar plates N/A=Not Applicable									
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	<i>10/22/12 HCL</i>	<i>11/2/12 HCL</i>	<i>11/2/12 HCL</i>						
<i>10/19/12 HCL</i>	<i>HAC</i>	<i>HEA</i>	<i>MSCD</i>				<i>10/22/12 HCL</i>	<i>N/A</i>	<i>N/A</i>
1	<i>NG</i>	<i>NG</i>	<i>NG</i>				<i>N/A</i>		
2									
3									
4									
5									
6									
7									
8									
9									
10									
11		<i>✓</i>	<i>✓</i>				<i>✓</i>		
12		<i>+</i>	<i>+</i>				<i>HEA</i> <i>MSCD</i> <i>> M</i>		
13		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
28 A		<i>+</i>	<i>+</i>				<i>HEA</i> <i>MSCD</i> <i>> M</i>		
28 B		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
28 C		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
28 D		<i>+</i>	<i>+</i>				<i>HEA</i> <i>MSCD</i> <i>> M</i>		
28 E	<i>✓</i>	<i>NG</i>	<i>NG</i>				<i>N/A</i>		
<i>N/A</i>									
<i>✓</i>									

ANALYST(S) *Hazlee B / Renew*

PAGE *22* OF *31* PAGES

Key: *+* = Pink colony *M* = Blackish gray mold-like colony *NG* = no growth

GENERAL CONTINUATION SHEET				PRODUCT <u>Methylpred. AC (PF) 80HG/HL</u>			SAMPLE NO. <u>764940</u>		
Results from Subculture of <u>+H10</u> Medium							N/A=Not Applicable		
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	<u>10/22/12 H10</u>	<u>11/2/12 H10</u>	<u>11/2/12 H10</u>						
<u>10/19/12 H10</u>	<u>HAC</u>	<u>HEA</u>	<u>MSCD</u>				<u>10/22/12 H10</u>	<u>N/A</u>	<u>N/A</u>
<u>29 A</u>	<u>NG</u>	<u>NG</u>	<u>NG</u>				<u>N/A</u>		
<u>29 B</u>		<u>NG</u>	<u>NG</u>				<u>N/A</u>		
<u>29 C</u>		<u>+</u>	<u>+</u>				<u>MEA</u> <u>MSCD</u> > <u>M</u>		
<u>29 D</u>		<u>NG</u>	<u>NG</u>				<u>N/A</u>		
<u>29 E</u>		<u>NG</u>	<u>NG</u>						
<u>30 A</u>		<u>NG</u>	<u>NG</u>						
<u>30 B</u>		<u>NG</u>	<u>NG</u>						
<u>30 C</u>		<u>+</u>	<u>+</u>				<u>MEA</u> <u>MSCD</u> > <u>M</u>		
<u>30 D</u>		<u>NG</u>	<u>NG</u>				<u>N/A</u>		
<u>30 E</u>		<u>NG</u>	<u>NG</u>						
<u>31 A</u>		<u>NG</u>	<u>NG</u>						
<u>31 B</u>		<u>NG</u>	<u>NG</u>						
<u>31 C</u>		<u>NG</u>	<u>NG</u>						
<u>31 D</u>		<u>NG</u>	<u>NG</u>						
<u>31 E</u>		<u>NG</u>	<u>NG</u>						
<u>32 A</u>		<u>NG</u>	<u>NG</u>						
<u>32 B</u>		<u>NG</u>	<u>NG</u>						
<u>32 C</u>		<u>NG</u>	<u>NG</u>						
<u>32 D</u>		<u>NG</u>	<u>NG</u>						
<u>32 E</u>		<u>NG</u>	<u>NG</u>						

ANALYST(S) Haydee B Romero

PAGE 23 OF 31 PAGES

Key: Y = Pink colony M = Blackish, gray mold-like colony

NG = No growth

GENERAL CONTINUATION SHEET	PRODUCT <u>Methy pred.AC(PF)80MG/ML</u> <u>Injectable</u>	SAMPLE NO. <u>764940</u>							
Results from Subculture of <u>TH10</u> Medium on agar plates N/A=Not Applicable									
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	10/22/12 <u>MAC</u>	11/2/12 <u>MEA</u>	11/2/12 <u>MSCD</u>						
33A	<u>NG</u>	<u>NG</u>	<u>NG</u>				<u>N/A</u>	<u>N/A</u>	<u>N/A</u>
33B		<u>NG</u>	<u>NG</u>				<u>N/A</u>		
33C		<u>+</u>	<u>+</u>				<u>MEA</u> <u>MSCD</u> > <u>M</u>		
33D		<u>NG</u>	<u>NG</u>				<u>N/A</u>		
33E									
34A									
34B									
34C									
34D									
34E									
35A									
35B									
35C									
35D									
35E									
36A									
36B									
36C									
36D									
36E	<u>✓</u>	<u>✓</u>	<u>✓</u>				<u>✓</u>		

ANALYST(S) Haydee B. Review

PAGE 24 OF 31 PAGES

Key: Y = Pink Colony M = Blackish, gray Mold-like Colony NG = No Growth

GENERAL CONTINUATION SHEET				PRODUCT <u>Methylpred. AC (PF) 80 MG/ML</u>				SAMPLE NO. <u>764940</u>	
Results from Subculture of <u>TH10</u> Medium									
on agar plates								N/A=Not Applicable	
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	10/22/12	11/2/12	11/2/12						
1/19/12 H2L	MAC	HEA	MSCD				10/22/12 H2L	N/A	N/A
37 A	NG	NG	NG				N/A		
37 B		NG	NG				N/A		
37 C		+	+				HEA MSCD > M		
37 D		+	+				HEA MSCD > M		
37 E		+	+				HEA MSCD > M		
38 A		NG	NG				N/A		
38 B			+						
38 C									
38 D									
38 E									
39 A									
39 B									
39 C									
39 D		✓	✓				✓		
39 E		+	+				HEA MSCD > M		
40 A		NG	NG				N/A		
40 B		+	+				HEA MSCD > M		
40 C		+	+				HEA MSCD > M		
40 D		NG	NG				N/A		
40 E		NG	NG				N/A		

ANALYST(S) Haydee B Roman

PAGE 2503 / PAGES

Key: Y = Pink colony

M = Blackish, gray mild-red colony

NG = No growth

ANALYST(S)

Haylee B Romero

PAGE 26 OF 31 PAGES

Key = Y = Pink Colony H = Blackish, gray Hotel-like colony

NG = No Growth

GENERAL CONTINUATION SHEET			PRODUCT <i>MeThy (pred. Ac (PF) 8046/HL Insectable</i>				SAMPLE NO. 764940		
Results from Subculture of <u>SCD</u> Medium on agar plates								N/A=Not Applicable	
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	10/22/12 <i>HL</i>	HEA	mscd						
10/19/12 <i>HL</i>	MAC	HEA	mscd				10/22/12 <i>HL</i>	N/A	N/A
1	NG	NG	NG				N/A		
2		NG	NG				N/A		
3		+	+				HEA mscd > Y		
4		+	+				HEA mscd > Y		
5		NG	NG				N/A		
6		+	+				HEA mscd > Y		
7		+	+				HEA mscd > Y		
8		+	+				HEA mscd > Y		
9		NG	NG				N/A		
10		+	+				HEA mscd > Y		
11		NG	NG				N/A		
12		+	+				HEA mscd > Y		
13		+	+				HEA mscd > Y		
28 A		+	+				HEA: Y mscd Y, M		
28 B		+	+				HEA mscd > M		
28 C		NG	NG				N/A		
28 D		+	+				HEA: Y mscd Y, M		
28 E	↓	+	+				HEA mscd > Y		
N/A									
↓									

ANALYST(S) *Haylee B. Rowen*

PAGE 27 OF 31 PAGES

Key: Y = pink colony M = blackish, gray mold-like colony NG = no growth

GENERAL CONTINUATION SHEET				PRODUCT <i>Methylpred. AC (PF) 8046/HL</i>			SAMPLE NO. <i>764940</i>		
Results from Subculture of <i>SCD</i> Medium on agar plates <i>N/A=Not Applicable</i>									
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	<i>10/22/12 HL</i> MAC	<i>4/2/12 HL</i> HEA	<i>4/2/12 HL</i> MSCD						
<i>10/19/12 HL</i>							<i>10/22/12 HL</i>	<i>N/A</i>	<i>N/A</i>
29 A	NG	+	+				HEA MSCD > Y		
29 B		NG	NG				N/A		
29 C		+	+				HEA: Y MSCD: Y, M		
29 D		+	+				HEA MSCD > Y		
29 E		NG	NG				N/A		
30 A		+	+				HEA MSCD > Y		
30 B		NG	NG				N/A		
30 C		+	+				HEA: Y MSCD: Y, M		
30 D		+	+				HEA MSCD > Y		
30 E		+	+				HEA MSCD > Y		
31 A		NG	NG				N/A		
31 B		+	+				HEA MSCD > Y		
31 C		+	+				HEA: Y, M MSCD: Y		
31 D		+	+				HEA MSCD > Y		
31 E		NG	NG				N/A		
32 A		+	+				HEA: Y, M MSCD: Y		
32 B		+	+				HEA MSCD > Y		
32 C		+	+				HEA: Y MSCD: Y, M		
32 D		+	+				HEA: Y MSCD: Y, M		
32 E	✓	NG	NG				N/A		

ANALYST(S) *Haylee B. Komer*PAGE *2* OF *3* PAGES

Key: Y = Pink colony

M = Blackish, gray mold-like colony

NG = No Growth

GENERAL CONTINUATION SHEET			PRODUCT <i>NeThylpred. AC(PF) 80MG/ML</i> <i>Injectable</i>			SAMPLE NO. <i>764940</i>			
Results from Subculture of <u><i>SCD</i></u> Medium on agar plates N/A=Not Applicable									
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	<i>10/28/12</i> MAC	<i>11/2/12</i> MEA	<i>11/2/12</i> MSCD						
<i>10/19/12</i> 33 A	<i>NG</i>	<i>+</i>	<i>+</i>				<i>10/22/12</i> HEA MSCD > Y	<i>✓</i> <i>IA</i>	<i>U/A</i>
33 B		<i>+</i>	<i>+</i>				HEA MSCD > Y		
33 C		<i>+</i>	<i>+</i>				HEA MSCD > Y		
33 D		<i>NG</i>	<i>NG</i>				N/A		
33 E		<i>NG</i>	<i>NG</i>				N/A		
34 A		<i>NG</i>	<i>NG</i>				N/A		
34 B		<i>NG</i>	<i>NG</i>				N/A		
34 C		<i>+</i>	<i>+</i>				HEA MSCD > Y		
34 D		<i>NG</i>	<i>NG</i>				N/A		
34 E		<i>NG</i>	<i>NG</i>				N/A		
35 A		<i>+</i>	<i>+</i>				HEA MSCD > Y		
35 B		<i>NG</i>	<i>NG</i>				N/A		
35 C		<i>+</i>	<i>+</i>				HEA : Y MSCD : Y & M		
35 D		<i>+</i>	<i>+</i>				HEA MSCD > Y, M		
35 E		<i>+</i>	<i>+</i>				HEA MSCD > Y		
36 A		<i>NG</i>	<i>NG</i>				N/A		
36 B		<i>+</i>	<i>+</i>				HEA : Y MSCD : Y & MC		
36 C		<i>NG</i>	<i>NG</i>				N/A		
36 D		<i>+</i>	<i>+</i>				HEA MSCD > Y		
36 E	<i>✓</i>	<i>NG</i>	<i>NG</i>				N/A		

ANALYST(S) *Haylee B. Romero*PAGE *29* OF *31* PAGES

Key: Y = Pink colony

M = Blackish, gray mild-like colony

NG = no growth

MC = Fine Cotton candy like, mild-like colony

Electronic Version of 431A

GENERAL CONTINUATION SHEET		PRODUCT <i>MeThylprod. AC (PF) 80MG/ML</i>				SAMPLE NO. 764940			
Results from Subculture of <u>SCD</u> Medium on agar plates N/A=Not Applicable									
SUB #	Agar used			Microscopic examination from			Colony Morphology	Subculture on	Biochemical identification
	10/22/12/10L	4/2/12/10L							
10/19/12 <i>HN</i>	<i>MAG</i>	<i>HEA</i>	<i>MSCD</i>				10/22/12/10L	<i>N/A</i>	<i>N/A</i>
37 A		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
37 B		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
37 C		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
37 D		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
37 E		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
38 A		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
38 B		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
38 C		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
38 D		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
38 E		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
39 A		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
39 B		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
39 C		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
39 D		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
39 E		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
40 A		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
40 B		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
40 C		<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		
40 D		<i>NG</i>	<i>NG</i>				<i>N/A</i>		
40 E	<i>✓</i>	<i>+</i>	<i>+</i>				<i>MEA</i> <i>MSCD</i> > <i>Y</i>		

ANALYST(S) *Haydee B. Pomeroy*
PAGE 30 **OF** 31 **PAGES**

Key: *Y* = Pink colony *M* = Blackish, gray mold-like colony

NG = No Growth

NRL.MSB.NECC

116 = no growth

Electronic Version of 431A

FDA E00277961

ATTACHMENT A

Sample # 764940

10/9/12 *WR*

From Direct inoculation of product into
MEA plate

36B





Photo taken by Jasna Brandt - Seawater 10/9/12



ATTACHMENT B

From direct inoculation of product onto
MEA plate

Sample# 764940
Date: 10/9/12
HR
37D



PHOTO TAKEN BY: Daniel J. Chalmers 10/9/12

Photo taken by Jesse Brant-Jacobsen 10/9/12

Sample# 764940
Date: 10/9/12
HR *mc*

37D

ATTACHMENT C
Photocopy of microscopic examination
of gray colony

Hagnibacter 100x



PHOTO TAKEN BY *Daniel J. Kelly* 10/9/12

ATTACHMENT D

Sample# 764940
Date: 10/9/12
HR *WZ*
37D

Photo/copy of microscopy examination of Pink Colony
Magnification 1000X

PHOTO TAKEN BY *Daniel J. Kelly* 10/9/12

Attachment E

Photocopy of microscope examination of pink colony

magnification 1000X

Sample# 764940

Date: 10/9/12

HR MC

37D

PHOTO TAKEN BY Danielle Dailly 10/9/12

Sample# 764940
Date: 10/9/12
HR
37D

ATTACHMENT F
From Direct Inoculation into MEA plate
magnification: 1000X



PHOTO TAKEN BY Daniel O'Leary 10/9/12


Product Unit Label

Sample # 764940
10/5/12 HR

Label removed from vial

METHYLPRED. AC (PF) 80MG/ML INJECTA
Lot# 08102012@51 Discard after 2/6/2013
SHAKE WELL ***SDV***
1 ML 764940 10/5/12 HR BCC

Label attached to foil bag




Product Unit Label

Sample # 764940
10/5/12 HR

Label removed from vial

METHYLPRED. AC (PF) 80MG/ML INJECTA
Lot# 08102012@51 Discard after 2/6/2013
SHAKE WELL ***SDV***
1 ML 764940 10/5/12 BCC

Label attached to foil bag




Product Unit Label

Sample # 764940
10/5/12 HR

Label removed from vial

METHYLPRED. AC (PF) 80MG/ML INJECTA
Lot# 08102012@51 Discard after 8/2013
SHAKE WELL***SDV***
1 ML 764940 10/5/12 HR

Label attached to foil bag



Food and Drug Administration Office of Regulatory Affairs

Collection Report

For Sample Number: 764940 13206132

This is an accurate reproduction of the original electronic record as of 12/10/2012

Flag	Flag Remarks				
Episode Number	Origin	Basis	Sample Type	FIS Smpl Num	Status
	Domestic	Other-Compln	Official	13206132	Completed
FEI	Date Collected	Product Code	Responsible Firm	PAC	Hours
3003623877	10/02/2012	64LDK45	Manufacturer	56R814	11
Compliance Num	Country of Origin				
	United States				
Related Smpl Num	Position Class	Sampling District	NDC Number	Permit Number	Storage Rqrmnt.
	INV	NOL-DO			Ambient
Dealer is Consumer	Crx/DEA Schedule	Recall Num	Consumer Compl. Num	Brand Name	
No				NECC	

Product Description

Methylprednisolone Acetate 80mg/ml packaged in 1ml clear glass vials. (Cont in Remarks)

Product Label

See continuation.

Reason for Collection

Collected during investigation of Aspergillus Meningitis outbreak.
Analyze for sterility (presence of Aspergillus) and chemical composition.

MFG Codes

"08102012@51"

Expiration Date

2/6/2013

Firm Legal Name	Address	Type of Firm	Firm FEI	FCE
St. Thomas Hospital	4220 Harding Pike Nashville, TN 37205-2005 US	Dealer	1000524946	
New England Compounding Center	697 Waverly St Framingham, MA 01702-8589 US	Manufacturer	3003623877	

Size of Lot	Est. Value	Rcpt Type	Carrier Name	Date Shipped
39 individual-1ml vials and 58-5 packs of 1ml vials		FDA484		

Description of Sample

Sample consists of 27 individual 1ml vials and 42-5 packs of vials.

Method of Collection

See continuation.

How Prepared

See continuation.

Collector's Identification on Package and/or Label

"764940 10/02/12 MHK" Subs 1-69

Collector's Identification on Seal

"764940 10/02/2012 Marshall H. Kinsey"

Sample Delivered To

UPS

Date Delivered

10/02/2012

Orig C/R & Records To

NOL-DO

Lab w/Split Sample

0

Lab

NRL

Document Number	Document Date	Document Type	Document Remarks
	10/02/2012	Other	Photographs of product labeling were taken.

Remarks

Date: 12/10/2012

Page: 1 of 3

**Food and Drug Administration Office of Regulatory Affairs
Collection Report**

For Sample Number: 764940

This is an accurate reproduction of the original electronic record as of 12/10/2012

See continuation.

Payment Amount	Payment Method	704(d) Sample	702(b) Portion	Collector's Name
\$0.00	No Charge	No	No	Marshall H Kinsey
Name of Signer	Date & Time of Signature			Meaning
Marshall H Kinsey	10/10/2012 10:24 AM ET			Collector

Food and Drug Administration Office of Regulatory Affairs

Collection Report

For Sample Number: 764940

This is an accurate reproduction of the original electronic record as of 12/10/2012

Continuation:

Product Label

1 ml vials are labeled in part, "****METHYLPRED. AC (PF) 80MG/ML INJECTA Lot#08102012@51***1 ML***necc****".

The foil packages are labeled in part, "NEW ENGLAND COMPOUNDING CTR***FRAMINGHAM, MA***Rx 1454300 JOHN W. CULCLASURE, MD 4230 HARDING RD. #901 NASHVILLE, TN***METHYLPRED. AC (PF) 80/MG INJECTABLE 5ML Lot#08102012@51***"

Bubble wrap bags are labeled in part, "****Total number of VIALS in this Bubble bag: *** Lot# 08102012@51".

Method of Collection

Using sterile gloves Inv. Kinsey reached into each of the four bubblewrap bags and individually removed the sample subs. Subs were placed in pre-labeled individual whirlpak bags which were held open by Inv. Smith. Whirlpaks were folded closed and placed in a clean cardboard box. Subs identified w/white self-adhesive labels & black print as shown. Subs 1-27 consist of individual vials collected from a single bubblewrap bag, subs 28-69 consist of 5-pack vials. Subs 28-41, 42-55, and 56-69 were collected from the remaining three bubblewrap bags.

How Prepared

Upon return to the office, subs were placed in a clean, unused, brown paper bag. The bag was folded closed, continuously wrapped in tape, and officially sealed as shown. Bag identified in black, handwritten ink as shown. Bag placed in a cardboard shipping container with bubblewrap and a completed FDA-525. Box prepared for shipment and taken to UPS by collector.

Remarks

(Cont from Product Description) Vials have a metal top with a white plastic lid. Vials are labeled with a white label w/ black and green print. Vials come packaged in metal foil packages, 5 per pack, which are labeled with white adhesive labels with black print. 5-pack vials come in bubblewrap bags, 20-5packs per bag. Bubblewrap bags are identified with white labels w/ black print.

A sealed FDA-525 containing a CD w/digital photos and the FDA-484 are attached to this collection report.

Approximately two thirds of products from available lots were collected. Dr. John W. Culclasure, Medical Director, stated that the hospital board had met and decided to retain 1/3 of all quarantined product for future testing by the hospital or other interested parties.

The dealer information listed in the Associated Firms section of this collection report is for the main hospital. This sample was collected from St. Thomas Outpatient Neurosurgical Center, 4230 Harding Rd. Plaza East Ste 901, Nashville, TN.